

Consensus Hip Systems



USA

Shalby Advanced Technologies, Inc.
1115 Windfield Way, Suite 100,
El Dorado Hills, CA 95762,
Ph: +1 916 355 7100

Singapore

Shalby Global Technologies Pte. Ltd.
Regus One Fullerton, 1 Fullerton Rd,
#02-01, One Fullerton, Singapore- 049231
Ph: +65 88707885 /+60 125065062

India

Mars Medical Devices Ltd.
Mondeal Heights, B-301 & 302,
Sarkhej - Gandhinagar Highway,
Ahmedabad, Gujarat-380015
Ph: +91 91671 20390



For over 20 years, the **Consensus Hip System** has provided proven fit and fill stems developed from well established design principles. Designed to help restore patient biomechanics.

The design philosophy of the Consensus system is to restore the original anatomic position of the femoral shaft to the acetabulum. The stem is designed to facilitate restoration of a wide variety of femoral geometries.

The Consensus stem incorporates a compound proximal wedge (CPW) design. The stem is tapered in the sagittal, coronal, and transverse plane. The effect of the CPW design improves load transfer for the cemented stem applications.

