



Taperset Hip System



USA

Shalby Advanced Technologies, Inc.
1115 Windfield Way, Suite 100,
El Dorado Hills, CA 95762,
Ph: +1 916 355 7100

Singapore

Shalby Global Technologies Pte. Ltd.
Regus One Fullerton, 1 Fullerton Rd,
#02-01, One Fullerton, Singapore- 049231
Ph: +65 88707885 /+60 125065062

India

Mars Medical Devices Ltd.
Mondeal Heights, B-301 & 302,
Sarkhej - Gandhinagar Highway,
Ahmedabad, Gujarat-380015
Ph: +91 91671 20390



The **TaperSet Hip System** features a reduced distal profile stem (RDP) which is designed to accommodate femurs with proximal-distal mismatch. The TaperSet RDP stem's reduced distal profile improves the proximal fit of the implant in the metaphysis. This unique reduced profile is the ideal choice to address a proximal/distal mismatch.

Taper-Set™

Reduced Distal Profile (RDP) Stems

