




# UniSyn Hip System



## USA

Shalby Advanced Technologies, Inc.  
1115 Windfield Way, Suite 100,  
El Dorado Hills, CA 95762,  
Ph: +1 916 355 7100

## Singapore

Shalby Global Technologies Pte. Ltd.  
Regus One Fullerton, 1 Fullerton Rd,  
#02-01, One Fullerton, Singapore- 049231  
Ph: +65 88707885 /+60 125065062

## India

Mars Medical Devices Ltd.  
Mondeal Heights, B-301 & 302,  
Sarkhej - Gandhinagar Highway,  
Ahmedabad, Gujarat-380015  
Ph: +91 91671 20390



[www.shalbytech.com](http://www.shalbytech.com)

[info@shalby.us](mailto:info@shalby.us)

Follow us on [f](#) [t](#) [in](#) [v](#)

[www.shalbytech.com](http://www.shalbytech.com)

The Consensus Orthopedics **UniSyn™ Hip System** has been designed with total flexibility in mind. All components are completely interchangeable. This allows for independent adjustment of both the size and orientation of the necks, bodies and stems.

**The UniSyn™ Hip System**  
**Completely interchangeable, limitless implant options**

All three components can be individually fit to specifically match the patients' anatomy. The UniSyn Hip was designed to address every surgical challenge encountered in primary or revision surgery.

Surgery is simplified by providing multiple surgical options and streamlined surgical techniques.



The UniSyn Hip System features patented<sup>9</sup> connection technology that allows the surgeon to independently fit the diaphysis and metaphysis, and then select the neck component that will best restore length, offset and anteversion of the patient's hip joint.



**Independent Neck and Body**

Allows for 360 degrees of version with the use of Modular Neck Vertical Lengths, or Lateral Offsets

**6 Proximal Body Options**

Multiple medial projection options for each cone size

**Multiple Stem Options**

Countless stem diameters and lengths to match any patients needs

**Distal flutes**

Improve rotational stability

**Coronal slot**

Reduces distal stem stiffness and helps reduce thigh pain