



Shalby
Advanced
Technologies, Inc.
Restoring Mobility, Improving Lives.

Revision Knee System

Instrument Catalog



Revision Knee System

Instrument Catalog

Table of Contents


Introduction.....	1
Tray 1	2
Tray 2	6
Tray 3	7
Tray 4	8
Tray 5	10
Tray 6 - Top.....	14
Tray 6 - Bottom.....	16
Tray 7	17
Tray 8	19

Introduction

Shalby Advanced Technologies, Inc Instrument sets are developed to provide surgeons with a simplified surgical experience. The Revision Knee System instrument set is a comprehensive system designed for improved ease-of-use and quick turnaround time for the entire OR staff. The Revision Knee instrument set allows a surgeon to revise any Consensus primary knee system.



Revision Knee System - Tray 1



Catalog #	Description	Quantity in set
0800-0-1600	Universal T-Handle (Quick Connector)	1



2800-0-0079	7.9mm, I.M. Drill Bit	1
-------------	-----------------------	---



2800-0-0132	3.2mm, Drill Bit	1
-------------	------------------	---



2800-0-0600	2 Piece Alignment Rods	1
-------------	------------------------	---



2800-0-1332	3in x 3.2mm, Fix. Pin with Cutout	4
2800-0-1432	4in x 3.2mm, Fix. Pin with Cutout	4



Revision Knee System - Tray 1

Catalog #	Description	Quantity in set
2801-0-0079	7.9mm, I.M. Drill Bit (Quick Connector)	1



2801-0-0132	3.2mm, Drill Bit	2
-------------	------------------	---




2810-0-0101	Short, I.M. Rod, Distal Femoral Alignment (Quick Connector)	1
-------------	--	---



2810-0-0350	1.27mm, Femoral Resection Indicator (Angle Wing)	1
-------------	--	---




2970-0-0027	.25in x .5in, Stop Drill Bit	2
-------------	------------------------------	---



GS 82.0240	G-source Pin Puller	1
------------	---------------------	---



Revision Knee System - Tray 1



Catalog #	Description	Quantity in set
2970-0-0007	Distal Femoral Cutting Guide	1



2970-0-0021	Femoral Alignment Guide	1
-------------	-------------------------	---



2970-0-1510	10mm x 150mm, Reamer Straight	1
2970-0-1512	12mm x 150mm, Reamer Straight	1
2970-0-1514	14mm x 150mm, Reamer Straight	1
2970-0-1516	16mm x 150mm, Reamer Straight	1
2970-0-1518	18mm x 150mm, Reamer Straight	1
2970-0-1520	20mm x 150mm, Reamer Straight	1
2970-0-1522	22mm x 150mm, Reamer Straight	1



2970-1-0008	I.M. Reamer Gauge	1
-------------	-------------------	---




Revision Knee System - Tray 1

Catalog #	Description	Quantity in set
2800-0-0500	Alignment Tower, Assembly	1



Revision Knee System - Tray 2



Catalog #	Description	Quantity in set
2981-0-0001	Left, Sz. 1, AP Cut Block	1
2981-0-0002	Left, Sz. 2, AP Cut Block	1
2981-0-0003	Left, Sz. 3, AP Cut Block	1
2981-0-0004	Left, Sz. 4, AP Cut Block	1
2981-0-0005	Left, Sz. 5, AP Cut Block	1
2981-0-0006	Left, Sz. 6, AP Cut Block	1



2975-0-0001	Left, Sz. 1, Femoral Trial	1
2975-0-0002	Left, Sz. 2, Femoral Trial	1
2975-0-0003	Left, Sz. 3, Femoral Trial	1
2975-0-0004	Left, Sz. 4, Femoral Trial	1
2975-0-0005	Left, Sz. 5, Femoral Trial	1
2975-0-0006	Left, Sz. 6, Femoral Trial	1



R11104	Femoral Impactor Extractor	1
--------	----------------------------	---

Revision Knee System - Tray 3



Catalog #	Description	Quantity in set
2982-0-0001	Right, Sz. 1, AP Cut Block	1
2982-0-0002	Right, Sz. 2, AP Cut Block	1
2982-0-0003	Right, Sz. 3, AP Cut Block	1
2982-0-0004	Right, Sz. 4, AP Cut Block	1
2982-0-0005	Right, Sz. 5, AP Cut Block	1
2982-0-0006	Right, Sz. 6, AP Cut Block	1




2976-0-0001	Right, Sz. 1, Femoral Trial	1
2976-0-0002	Right, Sz. 2, Femoral Trial	1
2976-0-0003	Right, Sz. 3, Femoral Trial	1
2976-0-0004	Right, Sz. 4, Femoral Trial	1
2976-0-0005	Right, Sz. 5, Femoral Trial	1
2976-0-0006	Right, Sz. 6, Femoral Trial	1



R11104	Femoral Impactor Extractor	1
--------	----------------------------	---



Revision Knee System - Tray 4



Catalog #	Description	Quantity in set
2971-0-0810	10mm x 80mm, Straight Stem Trial	1
2971-0-0812	12mm x 80mm, Straight Stem Trial	1
2971-0-0814	14mm x 80mm, Straight Stem Trial	1
2971-0-0816	16mm x 80mm, Straight Stem Trial	1
2971-0-0818	18mm x 80mm, Straight Stem Trial	1
2971-0-0820	20mm x 80mm, Straight Stem Trial	1
2971-0-0822	22mm x 80mm, Straight Stem Trial	1



2971-0-1110	10mm x 110mm, Straight Stem Trial	1
2971-0-1112	12mm x 110mm, Straight Stem Trial	1
2971-0-1114	14mm x 110mm, Straight Stem Trial	1
2971-0-1116	16mm x 110mm, Straight Stem Trial	1
2971-0-1118	18mm x 110mm, Straight Stem Trial	1
2971-0-1120	20mm x 110mm, Straight Stem Trial	1
2971-0-1122	22mm x 110mm, Straight Stem Trial	1

Revision Knee System - Tray 4


Catalog #	Description	Quantity in set
2971-0-1510	10mm x 150mm, Straight Stem Trial	1
2971-0-1512	12mm x 150mm, Straight Stem Trial	1
2971-0-1514	14mm x 150mm, Straight Stem Trial	1
2971-0-1516	16mm x 150mm, Straight Stem Trial	1
2971-0-1518	18mm x 150mm, Straight Stem Trial	1
2971-0-1520	20mm x 150mm, Straight Stem Trial	1
2971-0-1522	22mm x 150mm, Straight Stem Trial	1



2970-0-0003	Plug Taper Trial	1
-------------	------------------	---



Revision Knee System - Tray 5



Catalog # **Description** **Quantity in set**

2800-0-2110 Tibial Prep Guide Handle 1



2970-0-0011 5mm, Trial Screw Augment 6
2970-0-0012 10mm, Trial Screw Augment 4



2979-0-2001 Small, Peg Augment Trial 6
2979-0-2002 Medium, Peg Augment Trial 6
2979-0-2003 Large, Peg Augment Trial 4



2970-0-0010 0.25in Hex, Peg, Driver 1



2970-0-0026 3mm, Ball End Hex Driver 1



Revision Knee System - Tray 5


Catalog #	Description	Quantity in set
2800-0-0700	Femoral / Tibial Impactor Handle	1
2800-1-0701	Small Tibial / Femoral Head Impactor	1



2970-0-0022	Puller Assembly	2
-------------	-----------------	---



Revision Knee System - Tray 5



Catalog #	Description	Quantity in set
2983-0-0500	0 Lat, 5mm, Augment Tibial Trial	1
2983-0-0501	0 Med, 1 Lat, 5mm, Augment Tibial Trial	1
2983-0-0512	1 Med, 2 Lat, 5mm, Augment Tibial Trial	1
2983-0-0523	2 Med, 3 Lat, 5mm, Augment Tibial Trial	1
2983-0-0534	3 Med, 4 Lat, 5mm, Augment Tibial Trial	1
2983-0-0545	4 Med, 5 Lat, 5mm, Augment Tibial Trial	1
2983-0-0556	5 Med, 6 Lat, 5mm, Augment Tibial Trial	1
2983-0-0560	6 Med, 5mm, Augment Tibial Trial	1



2983-0-1000	0 Lat, 10mm, Augment Tibial Trial	1
2983-0-1001	0 Med, 1 Lat, 10mm, Augment Tibial Trial	1
2983-0-1012	1 Med, 2 Lat, 10mm, Augment Tibial Trial	1
2983-0-1023	2 Med, 3 Lat, 10mm, Augment Tibial Trial	1
2983-0-1034	3 Med, 4 Lat, 10mm, Augment Tibial Trial	1
2983-0-1045	4 Med, 5 Lat, 10mm, Augment Tibial Trial	1
2983-0-1056	5 Med, 6 Lat, 10mm, Augment Tibial Trial	1
2983-0-1060	6 Med, 10mm, Augment Tibial Trial	1




2979-0-0001	Sz. 1, Augment, Distal / Femoral Trial	2
2979-0-0002	Sz. 2, Augment, Distal / Femoral Trial	2
2979-0-0034	Sz. 3/4, Augment, Distal / Femoral Trial	2
2979-0-0056	Sz. 5/6, Augment, Distal / Femoral Trial	2

Revision Knee System - Tray 5

Catalog #	Description	Quantity in set
2979-0-1001	Sz. 1, Posterior, Augment Trial	2
2979-0-1002	Sz. 2, Posterior, Augment Trial	2
2979-0-1003	Sz. 3, Posterior, Augment Trial	2
2979-0-1004	Sz. 4, Posterior, Augment Trial	2
2979-0-1005	Sz. 5, Posterior, Augment Trial	2
2979-0-1006	Sz. 6, Posterior, Augment Trial	2



Revision Knee System - Tray 6 - Top



Catalog #	Description	Quantity in set
2970-0-0002	Body, Tibial Cut Block	1
2970-0-0009	Guide Reamer Proximal, Tibial	1
2970-0-0013	Bushing Guide, Tibial Cut Block	1
2970-0-0025	Top Tibial Alignment Guide	1



Revision Knee System - Tray 6 - Top

Catalog #	Description	Quantity in set
2980-0-0010	10mm, Bushing Guide	1
2980-0-0012	12mm, Bushing Guide	1
2980-0-0014	14mm, Bushing Guide	1
2980-0-0016	16mm, Bushing Guide	1
2980-0-0018	18mm, Bushing Guide	1
2980-0-0020	20mm, Bushing Guide	1
2980-0-0022	22mm, Bushing Guide	1



Revision Knee System - Tray 6 - Bottom




Catalog #	Description	Quantity in set
4410-20070-A	Femoral Head Impactor	1



2970-0-0020	Bottom Tibial Alignment Guide	1
-------------	-------------------------------	---




Revision Knee System - Tray 7



Catalog #	Description	Quantity in set
2943-1-1210	Left, Sz. 1-2, 10mm, PCL Trial Insert	1
2943-1-1212	Left, Sz. 1-2, 12mm, PCL Trial Insert	1
2943-1-1214	Left, Sz. 1-2, 14mm, PCL Trial Insert	1
2943-1-1216	Left, Sz. 1-2, 16mm, PCL Trial Insert	1
2943-1-1218	Left, Sz. 1-2, 18mm, PCL Trial Insert	1
2943-1-1220	Left, Sz. 1-2, 20mm, PCL Trial Insert	1
2943-1-1222	Left, Sz. 1-2, 22mm, PCL Trial Insert	1
2943-1-3410	Left, Sz. 3-4, 10mm, PCL Trial Insert	1
2943-1-3412	Left, Sz. 3-4, 12mm, PCL Trial Insert	1
2943-1-3414	Left, Sz. 3-4, 14mm, PCL Trial Insert	1
2943-1-3416	Left, Sz. 3-4, 16mm, PCL Trial Insert	1
2943-1-3418	Left, Sz. 3-4, 18mm, PCL Trial Insert	1
2943-1-3420	Left, Sz. 3-4, 20mm, PCL Trial Insert	1
2943-1-3422	Left, Sz. 3-4, 22mm, PCL Trial Insert	1
2943-1-5610	Left, Sz. 5-6, 10mm, PCL Trial Insert	1
2943-1-5612	Left, Sz. 5-6, 12mm, PCL Trial Insert	1
2943-1-5614	Left, Sz. 5-6, 14mm, PCL Trial Insert	1
2943-1-5616	Left, Sz. 5-6, 16mm, PCL Trial Insert	1
2943-1-5618	Left, Sz. 5-6, 18mm, PCL Trial Insert	1
2943-1-5620	Left, Sz. 5-6, 20mm, PCL Trial Insert	1
2943-1-5622	Left, Sz. 5-6, 22mm, PCL Trial Insert	1




Revision Knee System - Tray 7



Catalog #	Description	Quantity in set
2973-0-0001	Left, Sz. 1, Tibial Broach	1
2973-0-0002	Left, Sz. 2, Tibial Broach	1
2973-0-0003	Left, Sz. 3, Tibial Broach	1
2973-0-0004	Left, Sz. 4, Tibial Broach	1
2973-0-0005	Left, Sz. 5, Tibial Broach	1
2973-0-0006	Left, Sz. 6, Tibial Broach	1




2977-1-0001	Left, Sz. 1, Tibial Sizer	1
2977-1-0002	Left, Sz. 2, Tibial Sizer	1
2977-1-0003	Left, Sz. 3, Tibial Sizer	1
2977-1-0004	Left, Sz. 4, Tibial Sizer	1
2977-1-0005	Left, Sz. 5, Tibial Sizer	1
2977-1-0006	Left, Sz. 6, Tibial Sizer	1



Catalog #	Description	Quantity in set
2944-1-1210	Right, Sz. 1-2, 10mm, PCL Trial Insert	1
2944-1-1212	Right, Sz. 1-2, 12mm, PCL Trial Insert	1
2944-1-1214	Right, Sz. 1-2, 14mm, PCL Trial Insert	1
2944-1-1216	Right, Sz. 1-2, 16mm, PCL Trial Insert	1
2944-1-1218	Right, Sz. 1-2, 18mm, PCL Trial Insert	1
2944-1-1220	Right, Sz. 1-2, 20mm, PCL Trial Insert	1
2944-1-1222	Right, Sz. 1-2, 22mm, PCL Trial Insert	1
2944-1-3410	Right, Sz. 3-4, 10mm, PCL Trial Insert	1
2944-1-3412	Right, Sz. 3-4, 12mm, PCL Trial Insert	1
2944-1-3414	Right, Sz. 3-4, 14mm, PCL Trial Insert	1
2944-1-3416	Right, Sz. 3-4, 16mm, PCL Trial Insert	1
2944-1-3418	Right, Sz. 3-4, 18mm, PCL Trial Insert	1
2944-1-3420	Right, Sz. 3-4, 20mm, PCL Trial Insert	1
2944-1-3422	Right, Sz. 3-4, 22mm, PCL Trial Insert	1
2944-1-5610	Right, Sz. 5-6, 10mm, PCL Trial Insert	1
2944-1-5612	Right, Sz. 5-6, 12mm, PCL Trial Insert	1
2944-1-5614	Right, Sz. 5-6, 14mm, PCL Trial Insert	1
2944-1-5616	Right, Sz. 5-6, 16mm, PCL Trial Insert	1
2944-1-5618	Right, Sz. 5-6, 18mm, PCL Trial Insert	1
2944-1-5620	Right, Sz. 5-6, 20mm, PCL Trial Insert	1
2944-1-5622	Right, Sz. 5-6, 22mm, PCL Trial Insert	1



Revision Knee System - Tray 8



<i>Catalog #</i>	<i>Description</i>	<i>Quantity in set</i>
2974-0-0001	Right, Sz. 1, Tibial Broach	1
2974-0-0002	Right, Sz. 2, Tibial Broach	1
2974-0-0003	Right, Sz. 3, Tibial Broach	1
2974-0-0004	Right, Sz. 4, Tibial Broach	1
2974-0-0005	Right, Sz. 5, Tibial Broach	1
2974-0-0006	Right, Sz. 6, Tibial Broach	1



2978-1-0001	Right, Sz. 1, Tibial Sizer	1
2978-1-0002	Right, Sz. 2, Tibial Sizer	1
2978-1-0003	Right, Sz. 3, Tibial Sizer	1
2978-1-0004	Right, Sz. 4, Tibial Sizer	1
2978-1-0005	Right, Sz. 5, Tibial Sizer	1
2978-1-0006	Right, Sz. 6, Tibial Sizer	1





Shalby
Advanced
Technologies, Inc.
Restoring Mobility, Improving Lives.

USA

Shalby Advanced Technologies, Inc.
1115 Windfield Way, Suite 100,
El Dorado Hills, CA 95762,
Ph: +1 916 355 7100

Singapore

Shalby Global Technologies Pte. Ltd.
Regus One Fullerton, 1 Fullerton Rd,
#02-01, One Fullerton, Singapore- 049231
Ph: +65 88707885 /+60 125065062

India

Mars Medical Devices Ltd.
Mondeal Heights, B-301 & 302,
Sarkhej - Gandhinagar Highway,
Ahmedabad, Gujarat-380015
Ph: +91 91671 20390

